**Bedienungsanleitung
Instruction Manual
Instructions d'utilisation**

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Telefon: _____

Kaufdatum: _____

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CAUTION! To work with this microscope, sharp and pointed aids are being used. Please take care that this microscope and its accessories are stored at a place out of reach of children. Let children only work with this microscope under an adult's supervision! Keep packing material (plastic bags etc.) away from children!

Do not dispose of electric equipment together with household waste material!
In observance of European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE) and its implementation in accordance with national law, electric equipment that have reached the end of its life must be collected separately and returned to an environmentally compatible recycling facility. Discharged batteries and damaged re-chargeable batteries must be disposed of at special battery collection points. Information is available from your local disposal agent or local authority regarding the disposal of devices or batteries manufactured after the 01.06.2006

All parts

- | | |
|---------------------------------------|---|
| 1 dioptré scale | 12 object plate (glass) |
| 2 rubber eyecups | 13 blue filter disc for transmitted illumination |
| 3 eyepiece | 14 halogen lamp for transmitted illumination |
| 4 binocular viewing | 15 transmitted light dimmer |
| 5 microscope head | 16 incident light dimmer |
| 6 head holder | 17 cover of incident light source |
| 7 objective (2x / 4x) | 18 focusing wheel |
| 8 dust cover for objective | 19 power switch |
| 9 holding clips | 20 power supply (mains in) |
| 10 object plate (black/white plastic) | 21 fuse holder |
| 11 fixing screw for object plate | 22 halogen reflector lamp for incident illumination |
| | 23 blue filter disc for incident illumination |

1. Location

Before starting, choose a suitable location for microscoping. Please take care that your microscope stands on an even and solid ground. For observing, a 230 V / 50 Hz mains socket is required.

2. Illumination

With the Analyth ICD, you have three illumination options: incident light, transmitted light and both. Firstly, take care that the dimmers (15 and 16) are switched off (turn back in direction to the pointed corner of the triangle), while the power switch (19) is switched off. Connect the 230 V cord to the power supply (20) and then to the mains socket. Switch on the power switch (19). With the dimmer (15), you can adjust the brightness of transmitted light, and with the dimmer (16), you can adjust the brightness of incident light.

3. Observation

First, remove the dust cover (8) from the objective. There are two object plates for your choice: a black/white plastic plate (10, A) and a transparent glass plate (12, B). Insert the chosen object plate and fix it with the fixing screw (11). Place the object to be viewed in the centre of the object plate (10 or 12). Secure the object with clamps (9) if necessary.

For incident illumination, you can use both sides of the plastic object plate (10). If you have bright objects, the black side is advantageous; if you have dark objects, the white side is recommended for contrastive imaging.

For transmitted illumination, use the glass object plate (12). The blue filter disc (13) should not be removed in order to have the same colour temperature for incident and transmitted light, because incident illumination uses a blue filter disc (23), too. The objective's magnification factor (2x or 4x) can be adjusted by turning the objective shaft (7). It is recommended to start observations always with the lowest magnification (2x objective). In order to avoid incidence of disturbing light, put the rubber eyecups (2) upon both eyepieces (3). Now, look through the eyepieces (3). Adjust the sharpness of the image with the focusing wheel (18).

4. Adjustment of interpupillary distance

The Analyth ICD is equipped with an adjustable binocular viewing (4) on the microscope head (5). The individual interpupillary distance, or the distance between both eyepieces, is important for user's comfort. To set your interpupillary distance, move the eyepieces against each other by "folding" until you find the ideal position.

5. Dioptré adjustment

The compensation for the individual dioptrés of the eyes is done by the dioptré scale (1) on the right eyepiece. Adjust the dioptré scale on the right eyepiece to the middle position (silver circle). Shut your right eye and use the focusing wheel (18) to focus for your left eye. Now shut your left eye and adjust the focus for your right eye by turning the dioptré scale (1). Now, the microscope is adjusted for your eyes.

6. Exchanging the halogen reflector lamp (incident illumination)

Switch off the power switch (19) and pull out the plug of the power cord from mains socket. Wait until the hot lamp is cooled down certainly. For replacing the halogen reflector lamp, remove the cover of incident light source (17) (by turning it anti-clockwise) from the housing (G). Pull out straightly ahead the old halogen reflector lamp (22) from the lamp base (H). Insert a new halogen reflector lamp of the same type (12 V / 10 W). Reinstall the cover of incident light source (17) on the housing. Due to the danger of burning, never use the incident illumination without the cover (17 and 23)!

7. Exchanging the halogen lamp (transmitted illumination)

Switch off the power switch (19) and pull out the plug of the power cord from mains socket. Wait until the hot lamp is cooled down certainly. For replacing the halogen lamp, incline the microscope carefully to see the bottom side, loose the fixing screw of the lamp door (E) on the bottom side using a screw driver, and open the lamp door (F). Pull out straightly ahead the old halogen lamp (14) from the lamp base. Insert a new halogen lamp of the same type (12 V / 10 W). Reinstall the lamp door and tighten the fixing screw (E). Notice the properly touching of a new lamp; take care not to touch the lamp with bare fingers. E.g., use the protective envelope of the lamp, which should be delivered with it, or use a tissue, in order to grasp the bulb. This avoids contamination, which can reduce the brightness and life time of the lamp. If you have touched the new bulb with your bare fingers, clean it with alcohol before using it.

8. Exchanging the fuse

Switch off the power switch (19) and pull out the plugs of the power cord from mains socket and from mains in at the microscope (20). Open the fuse holder (reach with a small screwdriver into the loop and pull forwards). Take out the fuse holder (21). Exchange the old fuse by a new one of the same type (250 V / 2 A, F). Install the fuse holder und plug again.

9. Cleaning and maintenance

The microscope is a high quality optical appliance. Therefore you should make sure that dust or moisture do not come into contact with your microscope. Avoid putting fingerprints on any optical surfaces. If dust or dirt should nevertheless have gotten onto your microscope, or the accessories, you should first remove this with a soft brush. Then clean the soiled area with a soft, lint-free cloth. To remove fingerprints from the optical surfaces it is best to use a lint-free, soft cloth, to which you have applied some alcohol. After usage, you should stow the microscope and the accessories in their correct containers again in a dry and dust-free environment.

Remember: A well maintained microscope will keep its optical quality for years and thus maintain its value.

10. Technical data:

Binocular viewing

Eyepieces: 2 pieces WF 10x (diameter 30 mm)

Total magnification: 20x, 40x

(Total magnification is the product of the magnification factor of the objective and the magnification of the eyepiece.)

Incident illumination: halogen reflector lamp (12 V / 10 W)

Transmitted illumination: halogen lamp (12 V / 10 W)

11. Accessories

Dust protective cover

Lens paper

Fuse (250 V / 2 A, F)

Halogen reflector lamp (12 V / 10 W) - Halogen lamp (12 V / 10 W)

Object plate (glass)

Object plate (plastic black/white)

Blue filter disc

2 Eyecups

Eyepieces: 2 pieces WF 10x (diameter 30 mm)

Power cord

12. EEC conformity explanation

Meade Instruments Europe GmbH & Co. KG, resident in 46414 Rhede/Westf., Gutenbergstr. 2, Germany, explains the agreement with in the following specified EEC guidelines for this product:

EN 61326

EN 61000-3-2:2005

Product description: Stereo microscope
Microscope Model: Analyth ICD

Rhede, July 2, 2008

Meade Instruments Europe GmbH & Co. KG

Helmut Ebbert
Managing director

13. Warranty

The period of warranty is 2 years, beginning on the day of purchase. Please keep the cash receipt as evidence of purchase. Devices which become defective during the warranty period can be returned to the dealer where the device was bought. The repaired device or a new one will then be returned to you. In the case of defects which occur after the end of the warranty period, the devices can also be returned. However, repairs which become necessary after the end of the warranty period will be subject to a service fee.

Important:

Make sure to return the device carefully packed in the original packaging in order to prevent transport damage. Please also enclose the cash receipt (or a copy). This warranty does not imply any restriction of your statutory rights.

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