- DE** Gebrauchsanleitung
- GB** Instruction Manual
- FR** Manuel d'instruction
- ES** Manual de instrucciones

BRESSER®

Researcher ICD • Art. No. 58-03000 / 58-03100

VORSICHT!

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Kaufdatum:

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All parts (Fig. 1-5):

- ① Binocular head
- ② Objective holder with objective lenses
- ③ Eyepieces
- ④ Electric illumination
- ⑤ Microscope stage
- ⑥ Object clips
- ⑦ Focus wheel
- ⑧ Bulbs
- ⑨ Dioptre scale
- ⑩ Dust protective cover
- ⑪ Socket with fuse
- ⑫ Dimmer
- ⑬ Main switch
- ⑭ Source switch
- ⑮ Cover plate
- ⑯ Illumination housing
- ⑰ Rechargeable battery

Your Researcher ICD microscope comes with either Halogen (Art. No. 58-03000) or LED (Art. No. 58-03100) illumination. Operation of these models is slightly different.

1. Choosing a place for the microscope

Before starting, choose a sufficient location for microscoping. Please take care that your microscope stands on an even and solid ground. For observing, a 230 V/50 Hz wall connector is required. With LED illumination, a wall connector will be needed to recharge the internal batteries.

2. Illumination

With the Researcher ICD, you have three illumination options: upper light, transmission light and both. First, connect the 230 V cord to a wall connector. Now, engage the main switch (13). With the source switch (14), you can adjust the illumination type: Pos. I for upper light, pos. II for transmission light, pos. III for both. At position 0, both lamps are off. With the dimmer (12), one can adjust the brightness.

3. Observation

First, remove the dust protective cover (10) from the objective holder (2). Place the desired object in the middle of the stage (5) and secure it with the clamps (6) if necessary. The objective's magnification (2x or 4x) can be adjusted by turning the objective holder (2). It is recommended to start observations always with the lowest magnification. Now look through the eyepieces (3) at the binocular head (1). Adjust the sharpness with the focuser (7).

4. Focusing and dioptre compensation

Set the dioptre scale (9) to „0“ and observe an object. Now close your left eye and adjust the focus with the focuser (7) for your right eye. Then, close your right eye and adjust the focus for your left eye with the dioptre scale (9). Now, the microscope is adjusted for your eyes. Note the setting of the dioptre scale you need to use in future.

5. Eye distance

The Researcher ICD is equipped with an adjustable binocular head. To set your eye distance, move the eyepieces against each other until you found the ideal position.

6. Replacing the halogen bulb (Halogen version only)

First, disconnect the microscope from the 230 V source and allow the illumination to cool down (approx. 5 minutes). To change a defective Lamp (8), perform the following steps: For the bottom lamp, remove the covering glass (15). For the upper light, screw off the lamp housing (16) in order to pull out the lamp. Spare bulbs (12 V/10 W) can be purchased in electronic stores. After changing the bulb, clean the new bulb with alcohol before using it.

Note:

The illumination of the LED model has a life span of approx. 20,000 hours. In the case of a technical error please contact the service center.

7. Replacing the fuse (Halogen version only)

The Researcher ICD Halogen is equipped with two fuses (11) (250 V/0.5 A). First, disconnect the power cord. Screw off the fuse holder (11) with a fitting screwdriver. Replace the fuse by a new one (250 V/0.5 A).

8. Replacing the rechargeable batteries (LED version only)

In the case the batteries become defective, they can be replaced by new ones. Open the Battery section at the underside of the microscope and put three new rechargeable batteries (Type NiMH AA 1.2V, ca. 1200 mAh) in place (17). Please take care of the correct position!

CAUTION:

DO NOT use non-rechargeable batteries. This may cause fire and/or explosion!

9. Cleaning and maintenance

Your microscope is a high-quality optical instrument. Please take care that neither dust nor fingerprints take place on the optical surfaces. If necessary, clean the microscope body with a slightly moistened soft cloth. Use a separate clean cloth with some drops of alcohol to clean the optical surfaces. After using, the microscope should be protected with its dust cover.

Note:

Please consider that a sufficient treated microscope keeps its value and optical quality for many years.

10. Technical data

Binocular head

Eyepieces: 10x and 20x (2 pcs. each), 30.5 mm diameter

Objectives: 2x and 4x, Stereo

Total magnification: 20x, 40x, 80x (The total magnification can be calculated by multiplying the objective magnification with the eyepiece magnification.)

Illumination: 12 V/10 W or LED

11. EEC conformity explanation

Meade Instruments Europe GmbH & Co. KG, resident in 46414 Rhede/Westf., Gutenbergstr. 2, Germany, explains the agreement with in the following specified EEC guidelines for this product:

EN 55011:1999

EN 61547:1995

EN 61000-3-2:2000

EN 61000-3-3:1995

Product description: Incident-/Transmission-type microscope

Model: Researcher ICD 20x-80x

Rhede, 2007-06-01

Meade Instruments Europe GmbH & Co. KG

Helmut Ebbert

Managing Director

12. Only for EU countries

Do not dispose of electric equipment together with household waste material!

In observance of European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE) and its implementation in accordance with national law, electric equipment that have reached the end of its life must be collected separately and returned to an environmentally compatible recycling facility.

Discharged batteries and damaged re-chargeable batteries must be disposed of at special battery collection points. Information is available from your local disposal agent or local authority regarding the disposal of devices or batteries.

13. Warranty

The period of warranty is 2 years, beginning on the day of purchase. Please keep the cash receipt as evidence of purchase. Devices which become defective during the warranty period can be returned to the dealer where the device was bought. The repaired device or a new one will then be returned to you. In the case of defects which occur after the end of the warranty period, the devices can also be returned. However, repairs which become necessary after the end of the warranty period will be subject to a service fee.

Important:

Make sure to return the device carefully packed in the original packaging in order to prevent transport damage. Please also enclose the cash receipt (or a copy). This warranty does not imply any restriction of your statutory rights.

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© 2007 Meade Instruments Europe. No part of this manual may be reproduced, sent, transferred or be translated into another language in any form without written permission of the Meade Instruments Europe GmbH & Co. KG. Errors and technical changes excepted.