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2

(GB) A. Pump head mounting screw 2 off
 B. Rotating motor
 C. DC motor **RS** stock no. 255-9605 etc.
 D. Mounting mounting screws 4 off

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(F) C. Pince de tube

(I) C. Fermo del tubo

The pump head must be operated from the drive motor (RS stock no. 255-9605 etc.) by a rigid panel (up to 3mm thick) or the pre-drilled mounting plate, RS stock no. 330-828 with hole pattern (fig 1). Only the centre hole set is required if the motor and pump head are to be mounted on an existing panel.

Pump assembly

1. Mount the motor behind the panel using the four screws supplied (screw heads to be on the same side of the panel as the motor shaft).
2. Mount the pump head housing on the front side of the panel with the motor shaft collar locating through the circular hole in the housing base and the two threaded holes on the housing base coinciding with the two mounting holes on the panel. Now secure the housing with the two screws supplied (screws heads are on the motor side of the panel).
3. Turn motor shaft by hand until the flat section is adjacent to hole A on the side of the pump housing (Figure 3).
Now unscrew the hexagonal socket head grub screw on the roller assembly using key supplied and locate the assembly on the tip of the motor shaft so that the grub screw head is facing hole A.
Now positively push the roller assembly down onto the motor shaft ensuring that the grub screw remains aiming at the flat section of the motor shaft. It is useful to use a thin spacer e.g. a piece of metal at the base of the pump housing so that when the roller drops down along the motor shaft a gap of 0.5mm (approx.) is maintained between the roller and the housing for friction free rotation (fig 2.).
Tighten the grub screw through hole A.
Note: The two others screws on the roller assembly are preset to give optimum tube/roller performance and therefore should not be adjusted.

Silicone tube fitting (or replacement)

Note: Ensure that the tube inside the roller is evenly positioned half way up the depth of the housing.

Peristaltic pump operation

The roller occludes the tube which as its normal size draws in fluid which is then trapped by the two rollers and then expelled from the tube. It can be seen that this action is self priming and the pump can run dry.
The flow rate for a given tube internal diameter depends on the roller speed and therefore on the d.c. motor drive voltage. Without motor speed control it is impossible to achieve linear speed reduction easily by motor voltage reduction. Speed reduction beyond 50% however requires feedback control techniques to maintain the linearity of the voltage/speed characteristics.

Flow rate

This is specified for pumping water at 20°C with zero suction and delivery heads. 4.8mm internal diameter silicone tube (RS stock no. 330-834).

RS motor	Nom. speed (r.p.m.) at 12 Vd.c.	Flow rate ml/min. (1/hr.)
255-9598*	130	212 (12.72)
255-9605	65	106 (6.36)
255-9611	40	65 (3.9)

*This motor should only be used when higher flow rates are required, however, the operating life of the motor and the silicone tube are likely to be reduced to half the values obtained when using the 65 or 40 r.p.m. units. This is due to increased loading on the motor and the higher strain in the tubing.

Silicone tubing

RS stock no. 330-834 available in 3m length and can be cut to shorter lengths if required. Tube life depends on pumping speed, fluid temperatures, **suction and delivery pressures and the chemical compatibility between the tube and fluid (see table 1 below). Typical tube life when used with RS pump head with 65 r.p.m. motor is in the region of 300 hours at zero delivery and suction pressure. It is however recommended that the tube is regularly changed to avoid splitting and the fluid leaking into the pump.

**This should be in the range -20 to +100°C.

Table 1 - guide for fluids and solutions suitable for pumping through silicone tubing

Acetaldehyde	Melamine
Acetic acid, cold	Methyl alcohol
Acetone	Methyl isobutyl Ketone*
Alamask**	Milk
Ammonia gas	Mustard
Ammonium hydroxide	Natural gas
Aniline, cold	Nickel salts
Barium hydroxide	Nitrogen
Barium salts	Oils, animal, vegetable
Beer	Oxalic glycol
Beerwort	Oxygen, cold
Beet sugar liquors	Ozone
Blood	Penicillin
Boric acid	Phthalic acid
Calcium bisulphate	Potassium cyanide
Calcium hydroxide	Potassium dichromite
Calcium hypochlorite	Potassium salts
Calcium salts	Propyl alcohol
Cane sugar liquors	Rubber latex
Carbon dioxide	Silver nitrate
Carbon monoxide	Soap solutions
Carbonic acid	Soda ash
Caster oil	Sodium bicarbonate
Citric acid	Sodium bisulphate
Coconut oil	Sodium bisulphite
Cooper salts	Sodium borate
Corn oil	Sodium chloride
Detergent solutions	Sodium cyanide
Dextrone	Sodium dichromate
Dextrose	Sodium hydroxide
Diatomaceous slurry	Sodium hyochlorite
Dimethyl formamide	Sodium nitrite
Dioxane	Sodium silicate
Ethyl alcohol	Sodium sulphide
Ethylene	Sodium sulphite
Ethylene diamine	Sodium thiosulphate
Ethylene glycol	Soyabean oil
Ferric salts	Stannic chloride
Formaldehyde	Stannous chloride
Formic acid	Stearic acid
Gelatine	Sulphur
Glucose	Sulphur dioxide, dry
Glycerine	Sulphur trioxide, dry
Glycols	Tannic acid
Green sulphate liquor	Tartaric acid
Hydrogen peroxide	Isopropyl alcohol
Isobutyl alcohol	Tin chlorides
Kymeine	Urea
Lactic acid, cold	Vinegar
Lead acetate	Water, fresh, sea
Linseed oil	Whiskey
Magnesium salts	Wine
Mercury	Wort
Xylene*	Yeast

RS Components shall not be liable for any liability or loss of any nature (howsoever caused and whether or not due to RS Components' negligence) which may result from the use of any information provided in RS technical literature.

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